CareAssist[™] ES Hospital Bed

Simply practical





Simply practical

The CareAssist ES bed is our foundational hospital bed. We've designed it especially for busy caregivers who need to accomplish a multitude of tasks in a minimum amount of time. A number of one-and-done features distinguish the CareAssist ES bed as one of the most practical beds available.



Do more with easy, one-step features

- Bedside controls offer the **Dining Chair™** one-button feature
- Four-corner brake and steer for easy mobility
- Built-in bed extender lets a single caregiver add length to the bed
- Line-of-Site[™] angle indicator lets nurses monitor a patient's head angle at a glance



The centrally located **fifth wheel** adds maneuverability to let you steer the bed easily during transport.

Protect your patients from skin breakdown

Led by experienced clinical professionals, our comprehensive Safe Skin[™] program for the prevention and treatment of skin breakdown helps improve outcomes that are designed to be repeatable and sustainable.

The CareAssist ES bed contributes toward this effort with its SlideGuard[™] feature, which helps protect your patient's skin during repositioning by keeping his or her hips in place while the back of the bed elevates.

The bed's open architecture allows you to use a variety of surfaces, offering maximum flexibility to care for patients at every level of risk for skin breakdown.



Easy, one-step features

The CareAssist ES bed is designed to make bedside tasks easy and fast. Caregivers can enjoy the benefits of a bed that delivers essential and dependable technologies to meet the fundamental needs of the Acute Care setting.



One-button **Dining Chair** position enables caregivers to place patient in the chair position while the bed remains at a low height.



Point-of-Care[®] siderail controls are conveniently located in a central location, on both sides of the bed, where patient and caregiver interaction occurs.



Built-in bed extender allows caregivers to easily add length to the bed.



Dual-locking brake casters with an audible "brake not set" alarm that reminds caregivers to secure the bed.



With the **battery backup**, there's no need to find an outlet when repositioning the bed (high/low position, head-of-bed, knee section and foot angle) during patient transfer.



SlideGuard bed articulation technology helps support patient positioning.





The **Electric Vascular Foot** allows the caregiver to raise patient's feet without having to lift the weight of the patient's legs.



SafeView™ Alert System keeps staff informed of a patient's bed status.



The CareAssist ES bed connects to your hospital system to send data to **NaviCare™** Nurse Call, EMRs or other clinical applications.

Designed especially for busy caregivers who need to accomplish a **multitude of tasks** in a minimum amount of time



Line-of-Site Feature Built-in bed extender Open Architecture (Range of available mattresses) SlideGuard Technology Frame Ability to send data to Nurse Call, EMR or other Clinical Applications

Technical Specifications

I-Rom

an

•••

Frame width	102 cm (40")
Frame length with head/footboard (w/ bumpers)	228.60 cm (90″)
Sleep deck width normal length extended length	91.44 cm (36″) 203.20 cm (80″) 213.36 cm (84″)
Bed clearance under frame	19.05 cm (7.5″)
Optional scale system accuracy weigh in any position scale repeatability	1% of patient weight yes 0.2% of patient weight
Low position floor to deck	40 cm (15.75")
High position floor to deck	82.55 cm (32.5")
Max head elevation	65°
Max Trend and Reverse Trend	16°
Caster size	14.5 cm (5.7")
Weight limit (Safe working load)	226.8 kg (500 lbs)

Electric Vascular Foot

· SafeView System*

Brake Alarms

Hill-Rom is a leading global medical technology company with more than 6,500 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health. Hill-Rom people, programs, and product brands work towards one mission: Every day, around the world, we enhance outcomes for our patients and their caregivers.

mobility is life

Hill-Rom is committed to the concept that mobilising people early and supporting their independence, improves life - in the hospital and at home.

This product is not available in countries requiring a CE mark. This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments. These medical devices comply with all applicable regulations. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Medical beds manufacturer:

Hill-Rom SAS - B.P. 14 - Z.I du Talhouët - 56330 Pluvigner – France (also European Representative) Völker GmbH - Wullener Feld 79 - 58454 Witten - Germany Hill-Rom Holdings, Inc. - 1069 State Route 46 East - Batesville, IN 47006 - USA Lifts manufacturer: Liko AB - Nedre vägen 100 - 975 92 Luleå – Sweden

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products. ©2014 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. Doc. No: 5EN132302-01, 21-07-2014

Not all products/options are available in all countries. For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

